Foam Roller

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Fi

Eng: Compact Foam Roller

Used for self-massage of both muscles and tight structures of the body, also called "fascia". Often used by athletes as a complement or substitute for stretching and is also popular in rehabilitation. Begin rolling for a short time for each session. When you become more experienced, you can increase with both longer and more intervals.

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Eng: Back rolls

The exercise can free up tension in the thoracic spine which in turn often helps the neck and lower back to feel better.

Start as shown in the picture with the roll in the middle of the back. Place your hands as support under the neck and move the body in the direction towards the feet so that the roller is rolled up towards the neck. Stop before you reach the cervical vertebrae. Move the body so that the roller rolls downwards to the middle of the back. Try to relax and allow yourself to slightly stick out your chest during the performance. Slowly roll backwards and forwards for a few minutes. Spontaneous "cracks" that may arise from the back are usually not harmful. If the "crack" is consistent with sharp pain you should stop immediately and consult a physician!

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Eng: **Outside of thighs**

Suitable for pain on the outside of the knee. Also suitable for releasing tension in the outside of the thighs.

1. Place one thigh on the roller as shown in the picture. Start with the roller just above the knee.
2. Now move the body so that the roller rolls on the outside of the thigh up to the hip and back. Repeat for a few minutes. Switch legs.

Hint! You may feel some pain to begin with but as tensions release the discomfort will reduce.

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Eng: **Calf rolls**

Exercise for tired and tense calf muscles.

Place lower leg on the roller and the other foot next to it or the lower part of the leg as in the picture. Push yourself up a little with your hands so that your backside eases off the floor and now roll backwards and forwards at an easy pace from the Achilles tendon to the knee. Switch legs.

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Eng: Massages the body's soft tissues!

Compact Foam Roller is a cylinder, made of foam rubber with a specially designed pattern that provides deeper massage. In this way you can vary the intensity of the massage by using your own body weight.

- Patterned foam roller to stretch or massage your muscles.
- Increases the blood circulation.
- Improves mobility and alleviates muscle pain.
- Speeds up the recovery and reconstruction of the muscles.
- Free from latex and phthalates.

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